
Sonic Stim +

Ultrasound and TENS

Instructions For Use
Please read carefully

Dear Customer,

Thank you for choosing **Sonic Stim +**. TensCare stands for high-quality, thoroughly tested products for the applications in the areas of gentle electrotherapy, muscle toning, continence management and pain relief during labour.

Please read these instructions for use carefully and keep them for later use and observe the information they contain.

Best regards,

Your TensCare Team

Contents

1.0	INTRODUCTION	4
2.0	INTENDED USE	5
3.0	CONTRAINDICATIONS, WARNING & CAUTIONS	5
4.0	UNIT INFORMATION	10
5.0	INSTALLATION	11
5.1	BEFORE USE	11
5.2	CONNECTION	11
5.3	CONNECTION OF THE POWER ADAPTER	12
5.4	THERAPY MODE	12
5.5	DISCONNECTING	12
6.0	OPERATION	12
6.1.	MEASURES WITH REGARDS TO TREATMENTS	12
6.2.	OPERATING THE DEVICE	13
6.3.	THE APPLICATOR	17
7.0	MAINTENANCE	18
7.1.	CLEANING THE DEVICE	18
7.2.	CLEANING THE APPLICATOR	18
7.3.	CLEANING THE LEAD WIRE AND ADAPTER	18
7.4.	CLEANING THE ELECTRODE PAD	19
8.0	TROUBLESHOOTING	20
9.0	SPECIFICATIONS AND TECHNICAL DATA	21
10.0	STORAGE	22
11.0	EMC	23
12.0	DISPOSAL OF WASTE ELECTRICAL AND ELECTRONIC PRODUCTS (WEEE)	27
13.0	WARRANTY	28

Symbols Used

	Symbol for protection against electric shock: Type BF
	This symbol on the unit means "Refer to instructions for use".
	Temperature Limitation: indicates the temperature limits to which the medical device can be safely exposed.
	Serial Number: indicates the manufacturer's serial number so that a specific medical device can be identified.
	Disposal in accordance with Directive 2012/19/EU (WEEE)
	Catalogue Number: indicates the manufacturer's catalogue number so that the device can be identified.
	The name and the address of the manufacturer
	Date of manufacture
	Complies with the European Medical Device Directive (93/42/EEC) and amended by directive 2007/47/EC requirements. Notified body TUY Rheinland (CE01 97)
IPX7	Only for Ultrasonic head: Protected against the effects of temporary immersion in water
	Keep dry
	Class II symbol

1.0 Introduction

This manual has been written for the users of **Sonic Stim +**. It contains general information on the operation, precautionary practices, and maintenance information of the device. To maximize the use, efficiency, and life of the device, please read the manual thoroughly and become familiar with the controls, as well as the accessories before operating the device.

Pay attention to the following before using the **Sonic Stim +**:

- 1) Keep yourself informed of the contraindications (see chapter 3).
- 2) The device may not be used in proximity (i.e. less than 2 meters) to shortwave equipment.
- 3) The device may not be used in so-called "wet rooms" (hydrotherapy rooms).

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than those described in these operating instructions.

Therapy possibilities

Sonic Stim + is a therapy apparatus that offers both ultrasound therapy and electrotherapy. Pain affects the quality and enjoyment of life, especially for those who suffer chronic pain. **Sonic Stim +** is an ultrasound and electrotherapy device for the treatment of chronic and acute muscular pain. The applicator has a radiant surface of 4.0cm² and frequency of 1 MHz. Combination therapy of ultrasound and electrotherapy is ideal to localise trigger points and/or pain points.

The applicator

The ultrasound applicator for **Sonic Stim +** has one-frequency head. This applicator can now supply 1 MHz ultrasound. The head has excellent beam characteristics, meeting all requirements of the existing standards. The excellent beam characteristics, ergonomic design and effective contact control of the single-frequency applicator make optimal treatment possible.

2.0 Intended Use

Sonic Stim + is a Portable Ultrasound and TENS (Transcutaneous Electrical Nerve Stimulation) combination therapy device that generates deep ultrasonic waves within body tissues and TENS for the treatment of selected medical conditions such as symptomatic relief of chronic intractable pain, post-traumatic pain and post-surgical pain, muscle spasms, and joint contractures. **Sonic Stim +** is not recommended for the treatment of malignancies.

3.0 Contraindications, Warning & Cautions

The precautionary instructions found in this section and throughout this manual are indicated by specific symbols. Understand these symbols and their definitions before operating this equipment. The definitions of these symbols are as follows:

Caution: Text with a "CAUTION" indicator symbol will explain possible safety infractions that could have the potential to cause minor to moderate injury to an individual or damage to equipment.

Warning: Text with a "WARNING" indicator will explain possible safety infractions that will potentially cause serious injury to an individual and/or equipment damage.

Danger: Text with a "DANGER" indicator will explain possible safety infractions that are imminently hazardous situations that could result in death or serious injury.

Contraindications

- Do not use over or near bone growth centres (epiphyseal discs) until bone growth is complete.
- Do not use over a healing fracture.
- Do not use over the eyes.
- Do not use over the heart.
- Do not use over brain tissue.
- Do not use on patients with demand type cardiac pacemakers.
- Do not use on someone who is pregnant.
- Do not use on testicles.
- Do not use on patients post laminectomy.

- Do not use over an area of the spinal cord following a laminectomy.
- Do not use on areas of the body that lack sensation.
- Do not use on areas of post-traumatic sequelae.
- Do not use if the patient has an endoprosthesis / metal implants.
- Do not use on patients with implanted neurostimulation systems.
- Do not use to treat malignancies nor in the region where tumours or malignant tumours are present.
- Do not use on patients who have thrombophlebitis and/or varices.
- Do not use on patients experiencing septic inflammation.
- Do not use on patients who have diabetes mellitus.
- Do not use on patients who have osteoporosis.
- Do not use over ischemic tissues in patients with vascular disease where the blood supply would be unable to follow the increase in metabolic demand.
- Do not use over the carotid sinus nerves or arteries, laryngeal or pharyngeal muscles.
- Do not use on patients with hemorrhagic diatheses (excessive bleeding disorders).
- Do not use over areas that are under anaesthesia.
- Do not use on acute injuries
- Do not use on open wounds.
- Do not use if patient is feverish (pyrexia).
- Do not use on patient with tuberculosis.
- Do not use on patients who have localized inflammation.

Cautions

- Read, understand, and practice the precautionary and operating instructions. Know the limitations and hazards associated with using any ultrasound device. Observe the precautionary and operational decals placed on the unit.
- Keep informed of the contraindications.
- Do not operate the device when connected to any other medical device.
- Do not operate this device in an environment where other devices used intentionally radiate electromagnetic energy in an unshielded manner.
- Ultrasound should be routinely checked before each use to ensure that all controls function normally.
 - Check intensity control - make sure it properly adjusts the intensity of the ultrasonic power output in a stable manner.
 - Check treatment time control - make sure it terminates ultrasonic power output when the timer reaches zero.

- Do not use sharp objects such as pencil point or ballpoint pen to operate the buttons on the control panel.
- Handle the ultrasound applicator with care. Inappropriate handling of the Ultrasound applicator may adversely affect its characteristics.
- Before each use, inspect the Ultrasound Applicator for cracks to avoid the ingress of conductive fluid.
- Inspect Applicator cables and associated connectors before each use.
- The ultrasound therapy control unit is not designed to prevent the ingress of water or liquids. Ingress of water or liquids may cause malfunction of internal components of the device and therefore create risk of injury to the patient.
- Caution should be used:
 - with patients suspected or diagnosed with epilepsy.
 - with patients suspected or diagnosed with heart problems.
- Caution should be used in the presence of the following:
 - When there is a tendency to haemorrhage following acute trauma or fracture.
 - Following recent surgical procedures when muscle contraction may disrupt the healing process.
 - Over the menstruating or pregnant uterus.
 - Over areas of the skin which lack normal sensation.
- Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternate conductive medium or alternate electrode placement.
- Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.
- Never apply electrodes over irritated or broken skin.
- The device should be kept out of the reach of children.
- The device should be used only with the leads and electrodes recommended for use by the manufacturer.
- Do not use in the bath or shower. The device should not be submerged in water or other liquids as this will possibly damage the device and startle the patient.
- The use of heat and cold producing devices, such as electric heating blankets, heating pads or ice packs, may impair the performance of the electrodes or alter the patient's circulation and increase the risk of injury to the patient.
- The **Sonic Stim +** should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk for injury.

Warnings

- Care must be taken when operating this device around other equipment.
- Potential electromagnetic or other interference may occur to either this device or to the other equipment, or both. Minimise this interference by not using this device in conjunction with the other equipment.
- This device may not be used in proximity (i.e. less than 2 meters) to short-wave equipment.
- Avoid exposure to direct sunlight, rain, excessive dust, moisture, mechanical vibrations and shocks.
- This device may not be used in so-called "wet rooms" (hydrotherapy rooms).
- Before administering any treatment, you should become acquainted with the operating procedures for each program of treatment, as well as the indications, contraindications, warnings, and precautions. Consult other resources for additional information regarding the application of electrotherapy and ultrasound.
- Do not use solvents to clean this device.
- Do not use this device if it is damaged in any way.
- This device must only be serviced, repaired and opened by individuals at authorised sales centres.
- Dispose of this device in accordance with local regulations. Keep the operating instructions with the device.
- Pregnant and nursing women should use caution when using the device.
- Avoid use over or near bone growth centres until bone growth is complete.
- Treatment time should not exceed 30 minutes a day.
- Do not use a cell phone while operating the device.
- Patients with sensitivity to the coupling gel should use caution when using the device.
- Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
- Stimulation should not be applied over the neck, thorax, or mouth. Severe spasm of the laryngeal and pharyngeal muscles may occur, and the contractions may be strong enough to close the airway or cause difficulty in breathing.
- Stimulation should not be applied trans-cerebrally (across the head), over the carotid sinus (where the jaw meets the neck), over metal implants or in conjunction with sleep apnea or heart monitors.
- Stimulation should not be applied transthoracically. Since the introduction of electrical current into the heart may cause cardiac arrhythmias.
- Stimulation should not be applied to swollen, infected or inflamed areas or skin eruptions, e.g. phlebitis, thrombophlebitis, varicose veins, etc.

- Stimulation should not be applied over, or in proximity to, cancerous lesions.
- Always keep the ultrasound head in constant motion.
- Use ample conductive gel with the ultrasound head to ensure good coupling throughout the treatment. If needed, apply more when setting the intensity.
- Consult your doctor or physiotherapist if you are in any doubt whatsoever.

Danger

Patients with an implanted neurostimulation device must not be treated with or be in close proximity to any shortwave diathermy, microwave diathermy, therapeutic ultrasound diathermy, or laser diathermy anywhere on their body. Energy from diathermy (shortwave, microwave, ultrasound, and laser) can be transferred through the implanted neurostimulation system, can cause tissue damage, and can result in severe injury or death. Injury, damage, or death can occur during diathermy therapy even if the implanted neurostimulation system is turned "off."

Biohazardous Materials

Handle, clean, and dispose of components and accessories that have come in contact with bodily fluids according to national, local, and facility rules, regulations, and procedures.

Adverse Reaction


Skin irritation, inflammation, and electrode burns beneath the electrodes are potential adverse reactions. Perform the following procedures to avoid the negative effects of ultrasound therapy.

- **Applicator Movement** - If movement of the applicator is too slow, the patient may feel periosteal pain characterized by a deep ache or pain. If motion is too fast, or if the applicator does not maintain good contact with the skin, the therapeutic effect of the sound waves will be reduced, and the applicator may overheat.
- **Patient Susceptibility** - Some patients are more sensitive to ultrasound output and may experience a reaction similar to a heat rash. Be sure to inspect the treatment area during and following treatment and stop using if an adverse reaction occurs.

- **Coupling** - Coupling is described as contact between the applicator and the treatment site and may be accomplished through the use of a coupling agent, such as gel or lotion. Anything used as a coupling agent must be highly conductive. Air is a very poor conductor of ultrasonic waves.

4.0 Unit Information

- 1) TIME LED: indicates treatment times of **5** minutes, **10** minutes and **15** minutes
- 2) MODE LED: indicates intensity of the ultrasound Low (**L**), Medium (**M**) and High (**H**).
- 3) **TIME** button: Adjusts treatment times to either **5** minutes, **10** minutes and **15** minutes
- 4) **MODE** button: Adjusts the ultrasound intensity: Low (**L**), Medium (**M**) and High (**H**).
- 5) **+** button: Increases the intensity of stimulation
- 6) **-** button: Decreases the intensity of the stimulation
- 7) **STIM** LED: indicates the stimulation output state - when illuminate the stimulation is on.

- 8) POWER LED: indicates power state
- 9) **ON/OFF** switch: Power on by shifting up or power off by shifting down
- 10) Adapter connection point
- 11) Lead wire and electrode pad connection point
- 12) Ultrasound Head / Applicator
- 13) Main body

5.0 Installation

5.1 Before Use

Remove the device and all accessories from the box. Inspect the device for damages or missing parts and/or accessories. Report any damage or missing parts or accessories to your local dealer or retailer from which you purchased this device.

Kit Contents

- 1 x **Sonic Stim +**
- 1 x Instruction manual
- 1 x Electrode (50x100mm)
- 1 x Lead wire
- 1 x Power adapter
- 1 x Ultrasound transmission gel (85g)
- 1 x Quick start guide

5.2 Connection

Prior to connecting this device to the power supply, verify that the voltage and frequency stated on the rating label match the available power supply.

The power adapter is a part of the supply circuit on which the device's safety depends on. The approvals for the **Sonic Stim +** are only valid if used in combination with the type of adapter provided with the device.

Caution: It is not permitted to connect **Sonic Stim +** to any other type of adapter other than adapter provided.


Caution: Connection of accessories other than the ones specified by the manufacturer can adversely affect the safety of the patient and the proper function of the equipment; therefore, it is not permitted.

5.3 Connection of the power adapter

- 1) Connect the power adapter to the device's power cord.
- 2) Connect the power adapter to the wall outlet.

5.4 Therapy Mode

Sonic Stim + offers 2 treatment modes:

- 1) Combination: Ultrasound + Electrical Stimulation therapy
- 2) Ultrasound: Ultrasound therapy

5.5 Disconnecting

- 1) Turn off the **Sonic Stim +** by switching power **ON/OFF** switch from **ON** to **OFF** position.
- 2) Remove the power adapter from the wall outlet.

6.0 Operation

6.1. Measures with regards to treatments

Before treatment

- Ensure there are no contraindications to treatment.
- Clean the skin of the treatment area with soap or alcohol (70%).
- If the skin has excess hair, trim or shave hair for optimal treatment.
- Apply a reasonable amount of ultrasound transmission/conductive gel to the treatment area. Use only the ultrasound gel with CE mark.
- Ultrasonic Action Function Test: Place the probe horizontally, then apply several water drops on the surface of the probe, turn the device on and press the time button to activate the ultrasound device. You will be able to observe the ultrasonic action as the water droplets will appear to be dancing on the sound head and you may notice a small amount of steam being released. The water droplets on the probe start to perform one million vibrations per second showing the atomization phenomenon.

During treatment

- Move the ultrasound-head in a circular motion. The area treated should be twice the size of the diameter of the applicator.

- If experiencing poor transmission of ultrasound energy, it is advised to add more gel or reposition the ultrasound-head.

Caution: The ultrasound-head should be moved in a slow, flat, circular motion over the skin surface of the treatment area. Apply the sound head evenly (in time) over the treatment area; not too slow to avoid inducing heat; not too fast to prevent bad contact which would reduce the effectiveness of the treatment.

After treatment


- Clean the contact surface immediately after each treatment. Make sure that no ultrasound gel remains on the treatment head. We recommend cleaning the head and cable daily, using lukewarm water - do not immerse the device in water.
- The treatment heads can be disinfected using a cloth moistened with 70% alcohol.
- Check if there are any signs of improvement (e.g. pain, circulation or mobility).

6.2. Operating the device

6.2.1 Ultrasound Therapy

1) Apply Transmission Gel


Apply a layer of ultrasound transmission gel on the treatment area. The gel acts as a coupling substance and ensures effectiveness. The area treated should be two times the diameter of the treatment head.

Caution: Never apply the gel to the applicator. The applicator will register this as contact and may emit ultrasound energy, which could damage the applicator.


2) Switch on the device

Connect the power adapter according to section 5.3. Switch on the device, using power on/off switch by switching from **OFF** to **ON** position. The Power LED will light green.


3) Adjusting intensity

Press the **MODE** button to select the ultrasound intensity. The intensity has three levels, Low (**L**), Medium (**M**) and High (**H**), each level corresponds to a light indicator.


4) Adjusting treatment time

Press the **TIME** button to cycle through the treatment time (5, 10 and 15 minutes), as shown by the Time LED. When the time is chosen, the system will start working. During working time, the user can press the **TIME** button to adjust the treatment time.

5) Start treatment

Move the treatment head in a flat, slow, circular motion over the skin surface area that is covered with a layer of ultrasound transmission gel. Apply the head evenly (in time) over the treatment area.

Caution: The device has a load detection system for safety. When the treatment head does not have good contact with the skin,

the device will stop treatment automatically. During this time, the TIME LED will flash slowly (1 Hz). The device will not restore treatment until the contact is good.

Caution: The device has a temperature protection function. When the temperature of the treatment head exceeds 42°C the treatment will automatically stop and the TIME LED will flash quickly

(2Hz). The device will not restore treatment until the temperature is below 40°C.

6) Turn off the device



After the time duration has been completed, the device will automatically revert to standby state. Once your therapy session has been completed, turn off the device by sliding the **ON/OFF** switch downwards from **ON** to the **OFF** position.

6.2.2 Combination therapy

- 1) Connect the lead wire and electrode pad to the unit as shown by the illustrations to the right.


- Plug the lead wire into the connection point attached to the device,
- Connect the electrode pad with the lead wire.
- Make sure all connections are securely in place.

Caution: The device must be turned off before connecting the lead wires to the device.


- 2) Place electrode firmly on the skin after cleaning and drying the treatment area.

- Place the electrode pad on the area of the body indicated by your physician or therapist.
- Make sure the electrode pad is placed firmly to the skin and has made good contact between the skin and the pad.

- 3) Apply Transmission Gel


Apply a layer of ultrasound transmission gel on the treatment area. The gel acts as a coupling substance and ensures treatment effectiveness. The area treated should be two times the diameter of the treatment head

Caution: Never apply the gel to the applicator. The applicator will register this as contact and may emit ultrasound energy, which could damage the applicator.


- 4) Switch on the device

Connect the power adapter according to section 5.3. Switch on the device, using power **ON/OFF** switch by switching from **OFF** to **ON** position. The Power LED will light green.


5) Adjust ultrasound intensity

Press **MODE** button to select ultrasound intensity within low, middle and high. The LED will be lighted to indicate the selected intensity.


6) Adjust treatment time

Press the **TIME** button to cycle through the treatment time (5, 10 and 15 minutes), as shown by the Time LED indicators. When the time is chosen, the system will start working. During working time, the user can press the **TIME** button to adjust the treatment time.

7) Adjust stimulation intensity

Press the **+** button to increase the intensity of the stimulation. Press the **-** button to decrease the intensity of the stimulation. The **STIM** LED will flash every time the **+** or **-** button is pressed.

Note: The colour of the LED indicates the output stimulation intensity.

- Green light: Output intensity < 10 V;
- Orange light: Output intensity > 10 V.

Caution: The electrical stimulation cannot work without the ultrasound output. Users can only adjust the stimulation intensity after the ultrasound is working properly.

Caution: If the stimulation level is uncomfortable or becomes uncomfortable, reduce the stimulation intensity to a comfortable level and contact your medical practitioner if problems persist.

Caution: Move the ultrasound treatment head while you are adjusting the stimulation intensity to prevent the local skin temperature from becoming too high or burning.

Caution: Each step of increase output intensity is 1 V when the output intensity less than 5V; 0.5V/step when the output intensity is over 5V.


Caution: In the particular standard for Electrical Nerve and Muscle Stimulators, IEC 60601-2-10, it is recommended not to exceed a

current density of 2 mA r.m.s./cm², otherwise skin irritations or burns can occur. For current types that contain a DC component we recommend not to exceed a current density of 0.2 mA/cm².

8) Start treatment

Move the treatment head in a flat, slow, circular motion over the skin surface treatment area that is covered with a layer of ultrasound transmission gel. Apply the head evenly (in time) over the treatment area.

Caution: The device has a load detection system for safety. If the electrode pad or the ultrasound treatment head do not have

good contact with the skin, the STIM LED and TIME LED will flash and stop treatment after the output intensity of stimulation surpasses 5V. The intensity will automatically but slowly increase to the setting level after the pad and treatment head have made good contact with the skin.

Caution: The device has a temperature protection function. When the temperature of the treatment head exceeds 42°C the treatment will automatically stop and the TIME LED will flash quickly (2Hz). The device will not restore treatment until the temperature is below 40°C.

Caution: The device works without vibration. You must move the applicator with a slow but deliberate speed, flat against the treatment area and in a circular motion around the treatment area.

After finishing the treatment, the device will enter the waiting state. It is not recommended that the user restart treatment upon completion of therapy.

9) Turn off the device

After the time duration has been completed, the device will automatically revert to standby state. Once your therapy session has been completed, turn off the product by sliding the **ON/OFF** switch downwards from **ON** to the **OFF** position.

6.3. The applicator

The applicator is a precision instrument. Great care is taken in the development and production to obtain the best possible beam

characteristics. Rough treatment (jarring or dropping) can adversely affect these characteristics and must therefore be avoided.

7.0 Maintenance

7.1. Cleaning the device

Switch off the device and disconnect it from the power supply. The device can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution). If a more sterile cleaning is needed, use a cloth moistened with an antimicrobial cleaner.

Caution: Do not submerge the device in liquids. Should the unit accidentally become submerged, contact the dealer or Authorized Service centre immediately. Do not attempt to use the device that has been submerged in any liquid substrate until inspected and tested by a Service Technician Certified by an Authorized Service centre. Do not allow liquids to enter the ventilation holes.

7.2. Cleaning the applicator

The applicator should be regularly inspected for damage, e.g. hairline cracks, which could allow the penetration of liquids. Clean the contact surface immediately after each treatment. Make sure that no ultrasound gel remains on the applicator after using. We further recommend cleaning the head and cable daily, using lukewarm water. The applicator can be disinfected using a cloth moistened with 70% alcohol.

7.3. Cleaning the lead wire and adapter

Periodically wipe the lead wire and adapter clean with a cloth dampened with a mild soap solution, and then gently wipe them dry. Rubbing alcohol on the lead wire will damage the insulation and dramatically shorten their life.

7.4. Cleaning the electrode pad

- 1) Switch the power off and remove the pad from the skin and the lead wire.
- 2) Wash the pad when the adhesive surface becomes dirty and/or the pad is difficult to attach to the skin.
 - Wash the pad softly with your fingertips under slow running cold water for several seconds (do not use a sponge/cloth/sharp object like a nail on adhesive side, do not use detergents, chemicals or soap).
- 3) Dry the pads and let the adhesive surface air-dry completely (do not wipe with a tissue paper or cloth).
- 4) Replace the pad on the clear plastic film and store in plastic bag.

Caution: The life of the electrode pad may vary by the frequency of washing, skin condition, and storage state.

Caution: If the electrode pad no longer sticks to your skin or the electrode pad is broken, you should replace with a new electrode pad.

Caution: Before applying the electrode pad, it is recommended to wash and degrease the skin, and then completely dry the area.

Caution: Do not turn on the device when the electrode pad is not positioned on the body.

Caution: Never remove the electrode pad from the skin while the device is still powered on.

Caution: If a replacement electrode is necessary, use only electrode pads that are 2 inch x 4 inch (50*100mm), the same as the electrode pad provided with the **Sonic Stim +** device.

Caution: Use of electrode pads larger than provided may reduce the effect of the stimulation. Use of an electrode pad that is much smaller than the electrode pad provided with **Sonic Stim +** device may increase the chance of skin irritation or electrode burns occurring under the electrode pad.

Caution: Always use electrode pads that CE marked.

8.0 Troubleshooting

Note: If the following measures fail to alleviate the problem, please call your authorised agency or supplier.

Problem	Possible causes	Possible solution
LEDS fail to light up	Adapter contact failure	Ensure adapter is connected. Check the following contacts: All contacts are in place All contacts are not broken
Electrical stimulation is weak/cannot feel stimulation	Electrode pad has dried out or is contaminated	Replace with new electrode pad
	Electrode pad does not stick to skin well	Reconnect the electrode pad or replace with new electrode pad
	Lead wire old/worn/damaged	Replace with new lead wire
	Electrical stimulation intensity is low	Increase the output intensity
	Intensity is too high	Decrease intensity
Stimulation is uncomfortable	Electrode active area size is too small	Only 50x100mm electrode pads should be used
	Damaged or worn electrode or lead wire	Replace with new electrode pad or lead wire
	May not be operating the device according to the manual	Please check the manual before use
Stimulation stops	Poor electrode contact	Reapply electrodes, secure firmly
	Damaged or worn electrode or lead wire	Replace with new electrode pad or lead wire
	No contact medium	Use with appropriate ultrasound gel
Stimulation is ineffective	Improper electrode or applicator placement	Reposition electrode and/or applicator
	Unknown	Contact clinician

9.0 Specifications and Technical Data

Technical data of Ultrasound	
Frequency	1MHz \pm 10%
Ultrasound power control (MODE)	3 intensity levels(L, M, H)
Output power	3W(L), 4W(M), 5W(H)
Pulse rate	100Hz \pm 10%
Pulse width	10ms
Pulse duration	3ms-5ms
R_{tpa}	2-3.3
Duty factor	30%, 40%, 50%
Effective radiating area (AER)	4.0cm ²
Intensity	1.25w/cm ² (duty cycle at 50%)
R_{BN} (Max)	5.0
Beam type	Collimated
Waveform	Pulsed
Treatment time	5min, 10min, 15min
Material of treatment head	Aluminium

Technical data of Electrical Stimulation	
Treatment time	5min, 10min, 15min
Carrier Frequency (C.F.)	2.5kHz
Beat Frequency	1-120Hz
Output voltage	0-15V (500 Ω Load)
Stim power control	25 intensity levels

Technical data of Sonic Stim + device	
Service of life	2 years
Safety class	Class II, BF-Type
Dimensions	209 x 53 x 89mm (LxWxH)
Weight	235g

Technical data of power supply	
Supply voltage	AC 100-240V
Frequency	50/60Hz
Supply current	0.5A
Output voltage	DC 15V
Output current	1.2A
Dimensions	64 x 50 x 26.5mm (LxWxH)
Weight	120g

Environmental conditions		
Operating conditions	Temperature	5-40°C
	Relative humidity	30%-75%
	Atmospheric pressure	700-1060hPa
Storage and transport conditions	Temperature	-20-55°C
	Relative humidity	10%-93%
	Atmospheric pressure	700-1060hPa

10.0 Storage

For a prolonged pause in treatment, place the unit with the adapter, lead wire, electrode pad and manual back in the case. Store it in a dry room and protect it against heat, sunshine and moisture. Store the machine in a cool, well-ventilated place according to the storage conditions listed above. Never place any heavy objects on the machine.

11.0 EMC

The device needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information supplied in this manual.

Care must be taken when operating this device adjacent to or stacked with other equipment. Potential electromagnetic or other interference could occur to this or other equipment. Try to minimize this interference by not using other equipment in conjunction with it.

The performance of the device was determined to be essential performance. This device has been thoroughly tested and inspected to assure proper performance and operation.

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions		
The Sonic Stim + device is intended for use in the electromagnetic environment specified below. The customer or the user of the Sonic Stim + should ensure the device is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The Sonic Stim + device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Sonic Stim + device is suitable for use in all establishments including those directly connected to the public low-voltage power supply buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity			
The Sonic Stim + device is intended for use in the electromagnetic environment specified below. The customer or the user of the Sonic Stim + should ensure the device is used in such an environment.			
Immunity Test	IEC 6061 test level	Compliance	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±15kV air	±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line(s) to lines ±2kV line(s) to earth	±1kV line(s) to lines ±2kV line(s) to earth	Mains power quality should that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during mains power interruptions, the device should be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a

IEC 61000-4-8			commercial environment or hospital.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Table 3

Guidance and manufacturer's declaration – electromagnetic immunity			
The Sonic Stim + device is intended for use in the electromagnetic environment specified below. The customer or the user of the Sonic Stim + should ensure the device is used in such an environment.			
Immunity Test	IEC 6061 test level	Compliance	Electromagnetic environment guidance
Conducted RF IEC 610004-6	3 Vrms 150 kHz - 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any point of the Sonic Stim + device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d=1.2\sqrt{P}$</p> <p>$d=12\sqrt{P}$, 80MHz to 800MHz</p> <p>$d=2.3\sqrt{P}$, 800MHz to 2.7GHz</p> <p>where P is the maximum output power rating of the transmitter In watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p>
Radiated RF IEC 61000-4-3	10 V/m 80 MHz - 2.7 GHz	10 V/m	

			<p>Interference may occur in the vicinity of each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>Note: At 80 MHz - 800 MHz the higher frequency range applies.</p> <p>Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people</p>			
<p>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Sonic Stim + device is used exceeds the applicable RF compliance level above, should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Sonic Stim +.</p> <p>Over the frequency range 150 kHz - 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the Sonic Stim + device	
<p>The Sonic Stim + is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Sonic Stim + can help prevent electromagnetic interference by maintaining a minimum distance from portable and mobile RF communications equipment (transmitters) as recommended below, according to the maximum output power of the communications equipment.</p>	
Rated maximum	Separation distance according to frequency of transmitter / m

output power of transmitter / W	150 kHz – 80 MHz D=1.2√P	80 MHz – 800 MHz D=1.2√P	800 MHz – 2.7 GHz D=1.2√P
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz - 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

12.0 Disposal of Waste Electrical and Electronic Products (WEEE)

Fully discharged batteries must be disposed of in a specially labelled collection container, at toxic waste collection points or through an electrical retailer. You are under legal obligation to dispose of batteries correctly.

Please dispose of the device in accordance with the legal obligation in your area.

13.0 Warranty

Please contact your dealer or the device centre in case of a claim under the warranty. If you must send in the unit, enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **Sonic Stim +** is one year from date of purchase. In case of a warranty claim, the date of purchase must be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs completed under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded from the warranty:
 - All damage which has arisen due to improper treatment, e.g. no accordance with the user manual.
 - All damage which is due to repairing or tampering by customer or unauthorized third parties.
 - Damage which has arisen during transportation from the manufacturer to the consumer or to the service centre.

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TensCare aim to give you the best possible product and service. We listen to your suggestions and are constantly trying to improve our products. We also want to learn about the way our products are used, and the benefits they give. If you have anything you would like to share with us, please get in touch.

www.tenscare.co.uk

Follow us:

TensCare Ltd

TensCare Limited

@TensCareLtd

@tenscareltd

This manual is valid for the Sonic Stim + (CT1032)

This user manual is published by Shenzhen Dongdixin Technology Co., Ltd. Shenzhen Dongdixin Technology Co., Ltd. reserves the right to improve and amend it at any time without prior notice.

Amendments will however be published in a new edition of this manual.

All Rights Reserved.

CE
0197

EC Declaration of Conformity

Shenzhen Dongdixin Technology Co., LTD hereby declare that the products Sonic Stim + Portable Ultrasound and TENS Combo Device: CT1032 meet the provisions of Directive 93/42/EEC which apply to them. The medical device has been assigned class IIa according to Annex IX of the Directive 93/42/EEC. It bears the mark CE0197. The product concerned has been designed and manufactured under a quality management system according to Annex II, excluding section 4 of Directive 93/42/EEC. Compliance of the designated product with the Directive 93/42/EEC have been assessed and certified by the Notified Body: TÜV Rheinland LGA Products GmbH following the procedure relating to the EC Declaration of Conformity set out in Annex II, excluding section 4 of Company: Shenzhen Dongdixin Technology Co., LTD.

Distributed by:

TensCare Ltd.

9 Blenheim Road, Epsom,
Surrey, KT19 9BE, UK
+44 (0) 1372 723 434
info@tenscare.co.uk

Shenzhen Dongdixin Technology Co., Ltd.

Floor 1-2, No.3 Building, Fanshen Xusheng Industrial
Estate Xilixiaobaimang 518108 Nanshan District,
Shenzhen P. R. China

Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80, 20537 Hamburg Germany
Tel: 0049-40-2513175 Fax: 0049-40-255726