



SIRIUS XL LARYNGOSCOPE RANGE


EXCEPTIONAL QUALITY ● FIBRE OPTIC GLASS ● ENCLOSED CORES ● XENON & XLED

TIMESCO
INSTRUMENTS FOR LIFE



Highly polished enclosed, glass fibre core for focused illumination. The enclosed fibre optic cores, prevents any ingress of debris and allows easy cleaning.

The fibre core is positioned to offer the least obstructive view.



Atraumatic distal end.

Fully enclosed fibre optic core.



Clear product information etched on the blade.

3 metal ball bearings for secure fitting.

A green, fibre optic ISO 7376 fitting metal block.

Fully autoclaveable metal blade with an anti-glare satin finish. Covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.

EXCEPTIONAL QUALITY

The Sirius laryngoscope blades belong to the industry leading Sirius family of blades, exclusively available from Timesco. They are made from high grade stainless steel & incorporate an exceptional glass fibre optic core, allowing maximum light from Timesco Xenon or XLED laryngoscope handles.

FULLY AUTOCLAVABLE

The Sirius laryngoscope blades are fully autoclavable and available in Macintosh and Miller designs. These are the most popular designs and are only illustrated in this brochure.

FIVE YEAR GUARANTEE

The Sirius laryngoscope blades are supplied with an extensive 5 year guarantee; our mark of quality that you can trust.

PRODUCT SPECIFICATION



Fully autoclavable
[I.C.001/9](#)



CE marked



Latex free



ISO 7376 compatible



Supplied with multilingual instructions
[I.L010/6](#)




Compatible blades are colour coded



Extensive guarantee


MACINTOSH



Product Code	Description	Measurements
2955.150.05	Sirius XL Macintosh Size 0	83mm x 20mm
2955.150.10	Sirius XL Macintosh Size 1	95mm x 20mm
2955.150.15	Sirius XL Macintosh Size 2	115mm x 24mm
2955.150.20	Sirius XL Macintosh Size 3	135mm x 24mm
2955.150.25	Sirius XL Macintosh Size 4	160mm x 24mm
2955.150.30	Sirius XL Macintosh Size 5	176mm x 25mm

Macintosh designs feature a curved blade, are inserted into the vallecula, epiglottis fold, and moved laterally to expose the laryngeal cavity and vocal cords. These are the most popular blades worldwide.

MILLER



Product Code	Description	Measurements
2955.185.03	Sirius XL Miller Size 00	65mm x 18mm
2955.185.05	Sirius XL Miller Size 0	77mm x 18mm
2955.185.10	Sirius XL Miller Size 1	102mm x 18mm
2955.185.15	Sirius XL Miller Size 2	150mm x 19mm
2955.185.20	Sirius XL Miller Size 3	195mm x 19mm
2955.185.25	Sirius XL Miller Size 4	207mm x 22mm

Miller designs feature a straight blade, are inserted over the epiglottis and moved laterally to expose the laryngeal cavity and vocal cords. The smaller sizes 00, 0 and 1 are preferred by clinicians for paediatric and infant use.

In this brochure only our Sirius XL enclosed fibre optic light guide laryngoscope system is illustrated. We also manufacture the removable core fibre optic system, [Optima CLX](#), the standard light system Orion and the single use disposable system [Callisto](#), details on all are available on separate brochures.

Fully autoclaveable metal blade, covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.

The non-replaceable glass fibre bundle is fully enclosed in the blade which allows easy cleaning.

Rigid metal construction from fully autoclaveable, antimagnetic stainless steel.

Atraumatic distal end.

Tilting tip for elevation of the epiglottis, providing an improved view of the larynx, without excessive force, especially for difficult intubations.

Clear product information etched on the blade.

Lever positioned at the rear of the blade to control the tilting tip.

Green ISO 7376 fitting metal block compatible with fibre optic handles

3 metal ball bearings for secure fitting.



DIFFICULT INTUBATIONS


The Sirius XL Eclipse laryngoscope blades are an invaluable aid for difficult intubations as they feature a flexible tip for elevation of the epiglottis fold providing an improved view of the trachea without excessive force resulting in less patient trauma.



SIRIUS XL ECLIPSE BLADES

Description	Measurements
Sirius XL Eclipse Macintosh Size 2	115mm x 24mm
Sirius XL Eclipse Macintosh Size 3	135mm x 24mm
Sirius XL Eclipse Macintosh Size 4	160mm x 24mm

INTUBATION USING A TIMESCO SIRIUS XL ECLIPSE BLADE



The Timesco Sirius XL Eclipse laryngoscope blades benefit from an articulating tip to enable the elevation of the epiglottis fold; providing an improved view of the trachea without excessive force, resulting in less risk of patient trauma especially in difficult intubations. For example:

- Patients with restricted mouth opening.
- Patients with protruding teeth.
- Patients that have decreased neck movement.
- Patients with Spinal injuries.
- Patients with cervical neck collars in use.

Timesco also offer single use Eclipse blades in our Callisto range. Callisto Eclipse brochure: [L.L.069/1](#)

A polished metal collar reflects and directs light up towards the fibre core of the blade for maximum illumination.

Clear product branding with the product code and CE located on the opposite side to the Xenon logo. The LOT number is shown on the top of the handle.

Knurled handle design provides excellent grip.

Fully autoclaveable metal handle, covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.

Super-bright replaceable 2.5V Xenon bulb providing excellent, white illumination with a LUX output of 18,000. Bulbs provide up to 20 hours of use and a colour temperature of 2,900(k).

Colour coded band to clearly identify it is ISO 7376 fitting compatible with green fitting reusable and single use blades.

Powered by dry cell batteries or easily converted to a rechargeable handle by using the Timesco R Lithium Engines.

Secure screw cap for easy access to batteries.



XENON BULBS 2.5V

Product Code	Description
3000.900.20	Optima XL Bulbs Xenon 2.5V Box of 6
3000.900.21	Optima XI Bulb Xenon 2.5v Per

XENON HANDLES

Product Code	Description	
3000.350.05	Xenon Laryngoscope Handle, Small	146mm x 19mm
3000.350.10	Xenon Laryngoscope Handle, Medium	146mm x 29mm
3000.358.10	Xenon Laryngoscope Handle, Stubby	111mm x 32mm

A polished metal collar reflects and directs light up towards the fibre core of the blade for maximum illumination.

Clear product branding with the product code and CE located on the opposite side to the XLED logo. The LOT number is shown on the top of the handle.

Colour coded band to clearly identify it is ISO 7376 fitting compatible with green fitting reusable and single use blades.

Knurled handle design provides excellent grip.

Fully autoclavable metal handle (excluding lamp), covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.



RELIABLE ILLUMINATION...

Timesco's new and improved XLED laryngoscope handle is designed to outperform others and provide consistent, bright, reliable illumination when you need it most.

FLICKER FREE...

High performance and reliability from a clear, flicker free replaceable LED lamp. Intelligent electronics control the light intensity to be flicker free, even at low battery voltages.

BRIGHT LIGHT, EVERY TIME...

Improved specification gives a bright, controlled illumination at all times. The lux and colour temperature from the lamp is consistent throughout the lifespan of the lamp and the batteries.

NEW LAMP AT A GLANCE...

- Super-bright replaceable 3V LED lamp
- Flicker free, white illumination
- An exceptional lux output of 40,000 (at lamp)
- Long lasting and improved bulb life
- Colour temperature of 5,000(k)
- 3 x more lux output than other leading brands

Further details on the Optima XLED Handles are available in a separate Brochure: [B.L.040/1](#)



XLED BULBS

Product Code	Description
3000.950.20	XLED Laryngoscope Bulbs LED, Box of 6
3000.950.21	XLED Laryngoscope Bulb LED

XLED HANDLES

Product Code	Description	Measurements
3000.500.05	XLED Laryngoscope Handle, Small	146mm x 19mm
3000.500.10	XLED Laryngoscope Handle, Medium	146mm x 29mm
3000.508.10	XLED Laryngoscope Handle, Stubby	111mm x 32mm



TIMESCO ION-R RECHARGERS

Product Code	Description
3125.999.10	Timesco ION-R Charger (Single) Includes: Charger, power supply and both Timesco R inserts (AA and C)
3125.999.20	Timesco ION-R Charger (Dual) Includes: 2 x Chargers, 1 x power supply and 2 x each Timesco R inserts (AA and C)

Timesco ION-R Rechargers are multi voltage and are supplied with accessories for operation worldwide. For more details see our [ION-R brochure \(L.L.060/1\)](#).



Only the Timesco R Engine is inserted into the charger.

Both the small and medium Timesco R Engines can be charged in the single and double chargers.

Timesco R Engines can be fully charged in less than 2 hours.

The colour changing LED indicates the health of the battery and the charging cycle.

Side plates can be removed to daisy chain 3 chargers together to receive power from a single power supply.

Once Timesco R Engines are fully charged the recharger will switch to a low power mode to avoid using unnecessary electricity.

TIMESCO R BATTERIES

Product Code	Description
3115.150.20	Timesco R Engine - Small (AA Size)
3115.150.50	Timesco R Engine Medium (C Size)

Only one battery required per handle. Timesco supply small and medium rechargeable batteries to help you create a rechargeable system. For more details see our [Timesco-R brochure \(L.L.060/1\)](#).



Timesco produce Sirius XL Laryngoscope sets in hard cases. Listed are some of the popular sets, bespoke sets are also available on request.

Sirius XL Sets in Hard Durable Plastic Cases with Spare Bulb.

Hard cases supplied as appropriate to the contents of each set, example case shown only.

2955.975.80 Sirius XL Standard Set Fibre Optic 1 Medium XL Handle And 3 Blades Macintosh Nos. 2, 3 & 4 + Bulb In Case

Product Code	Description
2955.150.15	Sirius XL Macintosh Blade 2 Fibre
2955.150.20	Sirius XL Macintosh Blade 3 Fibre
2955.150.25	Sirius XL Macintosh Blade 4 Fibre
3000.350.10	Xenon Laryngoscope Handle Medium (C) Fibre
3000.900.21	Optima XI Bulbs Xenon 2.5v
2905.890.30	Plastic Laryngoscope Case For 3 Macintosh Blades + 1 Handle

2955.975.85 Sirius XL Standard Set Fibre Optic 1 Medium XLED Handle And 3 Blades Macintosh Nos. 2, 3 & 4 + Bulb In Case

Product Code	Description
2955.150.15	Sirius XL Macintosh Blade 2 Fibre
2955.150.20	Sirius XL Macintosh Blade 3 Fibre
2955.150.25	Sirius XL Macintosh Blade 4 Fibre
3000.500.10	XLED Laryngoscope Handle Medium (C) Fibre
3000.950.21	XLED Laryngoscope Bulb LED
2905.890.30	Plastic Laryngoscope Case For 3 Macintosh Blades + 1 Handle

2955.985.70 Sirius XL Paediatric Set 1 Small XL Handle and 2 Blades Miller Nos. 0 & 1 + Bulb In Case

Product Code	Description
2955.185.05	Sirius XL Miller Blade 0 Fibre
2955.185.10	Sirius XL Miller Blade 1 Fibre
3000.350.05	Xenon Laryngoscope Handle Small (AA) Fibre
3000.900.21	Optima XI Bulbs Xenon 2.5v
2905.890.42	Plastic Laryngoscope Case For 2 Miller Blades + 1 Handle

2955.989.95 Sirius XL Combination set 2 Handles 1 x Medium XL, 1 x Small XL and 5 Blades Macintosh Nos. 2, 3 & 4, Miller Nos. 0 & 1 + Bulb In Case

Product Code	Description
2955.150.15	Sirius XL Macintosh Blade 2 Fibre
2955.150.20	Sirius XL Macintosh Blade 3 Fibre
2955.150.25	Sirius XL Macintosh Blade 4 Fibre
2955.185.05	Sirius XL Miller Blade 0 Fibre
2955.185.10	Sirius XL Miller Blade 1 Fibre
3000.350.05	Xenon Laryngoscope Handle Small (AA) Fibre
3000.350.10	Xenon Laryngoscope Handle Medium (C) Fibre
3000.900.21	Optima XL Bulbs Xenon 2.5v
2905.890.75	Plastic Laryngoscope Case For 3 Mac & 2 Mill Blades + 2 Handles



SIRIUS XL ECLIPSE - INDIVIDUAL BLADES

Product Code	Description	Measurements
2955.175.15	Sirius XL Eclipse Macintosh 2 in case	115mm x 24mm
2955.175.15X	Sirius XL Eclipse Macintosh 2 without case	115mm x 24mm
2955.175.20	Sirius XL Eclipse Macintosh 3 in case	135mm x 24mm
2955.175.20X	Sirius XL Eclipse Macintosh 3 without case	135mm x 24mm
2955.175.25	Sirius XL Eclipse Macintosh 4 in case	160mm x 24mm
2955.175.25X	Sirius XL Eclipse Macintosh 4 without case	160mm x 24mm



SIRIUS XL ECLIPSE BLADES

Product Code	Description
2955.177.20	Sirius XL Eclipse Macintosh 2 and 3 in case



SIRIUS XL ECLIPSE BLADES AND XENON HANDLE

Product Code	Description	Measurements
2955.176.15	Sirius XL Eclipse Macintosh 2 and Xenon Handle in case	115mm x 24mm
2955.176.20	Sirius XL Eclipse Macintosh 3 and Xenon Handle in case	135mm x 24mm
2955.176.25	Sirius XL Eclipse Macintosh 4 and Xenon Handle in case	160mm x 24mm



SIRIUS XL ECLIPSE BLADES AND XENON HANDLE

Product Code	Description
2955.177.25	Sirius XL Eclipse Macintosh 3 & 4 blades & Xenon Handle set in case



SIRIUS XL ECLIPSE BLADES AND XLED HANDLE

Product Code	Description	Measurements
2955.178.15	Sirius XL Eclipse Macintosh 2 Blade & XLED Handle in case	115mm x 24mm
2955.178.20	Sirius XL Eclipse Macintosh 3 Blade & XLED Handle in case	135mm x 24mm
2955.178.25	Sirius XL Eclipse Macintosh 4 Blade & XLED Handle in case	160mm x 24mm



SIRIUS XL ECLIPSE BLADES AND XLED HANDLE

Product Code	Description
2955.179.25	Sirius XL Eclipse Macintosh 3 & 4 Blades & XLED Handle Set in case



SIRIUS XL ECLIPSE BLADES AND XENON HANDLE

Product Code	Description
2955.180.50	Sirius XL Eclipse Macintosh 2, 3 & 4 blades & XENON Handle Set in case



SIRIUS XL ECLIPSE BLADES AND XLED HANDLE

Product Code	Description
2955.182.50	Sirius XL Eclipse Macintosh 2, 3 & 4 blades & XLED Handle Set in case



SIRIUS XL ECLIPSE - EMPTY CASE (2 BLADES)

Product Code	Description
2999.175.75	Plastic Case for 2 Sirius XL Eclipse Blades and Handle



SIRIUS XL ECLIPSE - EMPTY CASE (3 BLADES)

Product Code	Description
2999.185.50	Plastic Case for 3 Sirius XL Eclipse Blades and Handle

Product Code	Description	Dimensions
SIRIUS MACINTOSH LARYNGOSCOPE BLADES		
2955.150.05	Sirius XL Macintosh Size 0	83mm x 20mm
2955.150.10	Sirius XL Macintosh Size 1	95mm x 20mm
2955.150.15	Sirius XL Macintosh Size 2	115mm x 24mm
2955.150.20	Sirius XL Macintosh Size 3	135mm x 24mm
2955.150.25	Sirius XL Macintosh Size 4	160mm x 24mm
2955.150.30	Sirius XL Macintosh Size 5	176mm x 25mm
SIRIUS MILLER LARYNGOSCOPE BLADES		
2955.185.03	Sirius XL Miller Size 00	65mm x 18mm
2955.185.05	Sirius XL Miller Size 0	77mm x 18mm
2955.185.10	Sirius XL Miller Size 1	102mm x 18mm
2955.185.15	Sirius XL Miller Size 2	150mm x 19mm
2955.185.20	Sirius XL Miller Size 3	195mm x 19mm
2955.185.25	Sirius XL Miller Size 4	207mm x 22mm
XENON LARYNGOSCOPE HANDLES		
3000.350.05	Xenon Laryngoscope Handle, Small	146mm x 19mm
3000.350.10	Xenon Laryngoscope Handle, Medium	146mm x 29mm
3000.358.10	Xenon Laryngoscope Handle, Stubby	111mm x 32mm
XLED LARYNGOSCOPE HANDLES		
3000.500.05	XLED Laryngoscope Handle, Small	146mm x 19mm
3000.500.10	XLED Laryngoscope Handle, Medium	146mm x 29mm
3000.508.10	XLED Laryngoscope Handle, Stubby	111mm x 32mm
OPTIMA XL BULBS XENON		
3000.900.20	Optima XL Bulbs Xenon 2.5V	Box of 6
3000.900.21	Optima XI Bulb Xenon 2.5v Per	
OPTIMA XLED BULBS LED		
3000.950.20	XLED Laryngoscope Bulbs LED	Box of 6
3000.950.21	XLED Laryngoscope Bulb LED	
TIMESCO BATTERIES		
3115.150.20	Timesco R Engine Small (AA Size)	Small (AA Size)
3115.150.50	Timesco R Engine Medium (C Size)	Medium (C Size)
TIMESCO ION-R CHARGER		
3125.999.10	Timesco ION-R Charger (Single) Includes: Charger, power supply and both Timesco R inserts (AA and C)	
3125.999.20	Timesco ION-R Charger (Dual) Includes: 2 x Chargers, 1 x power supply and 2 x each Timesco R inserts	

LARYNGOSCOPE SETS

Product Code	Description	Measurements
SIRIUS XL ECLIPSE BLADES IN CASE		
2955.175.15	Sirius XL Eclipse Macintosh 2 in case	115mm x 24mm
2955.175.15X	Sirius XL Eclipse Macintosh 2 without case	115mm x 24mm
2955.175.20	Sirius XL Eclipse Macintosh 3 in case	135mm x 24mm
2955.175.20X	Sirius XL Eclipse Macintosh 3 without case	135mm x 24mm
2955.175.25	Sirius XL Eclipse Macintosh 4 in case	160mm x 24mm
2955.175.25X	Sirius XL Eclipse Macintosh 4 without case	160mm x 24mm

SIRIUS XL ECLIPSE BLADES AND HANDLE IN CASE

2955.176.15	Sirius Eclipse Macintosh 2 and Xenon Handle in case	115mm x 24mm
2955.176.20	Sirius Eclipse Macintosh 3 and Xenon Handle in case	135mm x 24mm
2955.176.25	Sirius Eclipse Macintosh 4 and Xenon Handle in case	160mm x 24mm
2955.178.15	Sirius Eclipse Macintosh 2 Blade & XLED Handle in case	115mm x 24mm
2955.178.20	Sirius Eclipse Macintosh 3 Blade & XLED Handle in case	135mm x 24mm
2955.178.25	Sirius Eclipse Macintosh 4 Blade & XLED Handle in case	160mm x 24mm
2955.177.20	Sirius XL Eclipse Macintosh 2 and 3 in case	
2955.177.25	Sirius Eclipse Macintosh 3 & 4 Blades & Xenon Handle set in case	
2955.179.25	Sirius Eclipse Macintosh 3 & 4 Blades & XLED Handle Set in case	
2955.180.50	Sirius Eclipse Macintosh 2, 3 & 4 Blades & Xenon Handle Set in case	
2955.182.50	Sirius Eclipse Macintosh 2, 3 & 4 Blades & XLED Handle Set in case	

SIRIUS XL ECLIPSE - EMPTY CASE

2999.175.75	Plastic Case for 2 Eclipse Blades and Handle	
2999.185.50	Plastic Case for 3 Eclipse Blades and Handle	

SIRIUS XL LARYNGOSCOPE SET IN CASE

2955.975.80	Sirius XL Standard Set Fibre Optic 1 Medium XL Handle and 3 Blades Macintosh Nos. 2, 3 & 4 + Bulb In Case	
2955.975.85	Sirius XL Standard Set Fibre Optic 1 Medium XLED Handle and 3 Blades Macintosh Nos. 2, 3 & 4 +Bulb In Case	
2955.985.70	Sirius XL Paediatric Set 1 Small XL Handle and 2 Blades Miller Nos. 0 & 1 + Bulb in Case	
2955.989.95	Sirius XL Combination set 2 Handles 1 x Medium XL Handle, 1 x Small XL Handle and 5 Blades Macintosh Nos. 2, 3 & 4, Miller Nos. 0 & 1 + Bulb In Case	

Contact Us

find out more today

For over 55 years' Timesco have been at the forefront of Medical product design and innovation. As a manufacturer and supplier of laryngoscope products, we have an unrivalled portfolio supplying both the UK directly and internationally through distributors.

As a family owned business, with headquarters in Basildon England, we strive to deliver products made to an exceptionally high standard; exceeding internationally recognised quality standards and customer expectations. We are recognised as one the most progressive and innovative companies within the industry and have a dedicated product development team which focuses on the changing needs of our clients and advances in clinical practice.

Contact Timesco today for further information or to request product samples.

t +44 (0) 1268 297710

e sales@timesco.com

W timesco.com



ISO 13485



CE 1639



FDA Manufacturers
License (USA)

(Registration Number: 8021788)



SFDA
CHINA



Ministry of Health
Ministry of Health

Classification: TIMESCO products are classified under the Harmonised Tariff: 90-1890-8400.

Timesco Healthcare Ltd, Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | **E** info@timesco.com | timesco.com

Ref:B.L.045/1 E&O.E. Timesco Healthcare Ltd. 11/20

TIMESCO
INSTRUMENTS FOR LIFE