



# Enhance 270 Hoist

Effortless transfers, maximum confidence

Providing  
comfortable,  
secure transfers  
for users weighing  
up to 272kg



**drive**

**DeVilbiss**  
HEALTHCARE

[drivedevilbiss.co.uk](http://drivedevilbiss.co.uk)

# Enhance 270 Hoist

The Enhance 270 Hoist is a smooth-operating, electrically powered mobile patient lift, designed to support and promote safe patient handling and transfer for both the user and caregiver.

Effortlessly manoeuvrable and stable thanks to its fully powered operating base, the Enhance 270 provides comfortable and secure transfers for users weighing up to 272kg (43st).



Boom

Six point cradle offers wide footprint

Emergency manual lowering mechanism

Emergency stop button

Ergonomic push handles for comfortable manoeuvring

Battery/control pack can be charged on or off board

Electrically adjustable leg operation

Easy roll front castors for optimal positioning

Large, braked rear castors provide good control and stability

A Drive DeVilbiss Smart/Safety control module is available to purchase separately. This stores useful servicing information about the lift that can be recalled when required.

- Accurate service data available at the touch of a button
- Optimised product lifetime
- Ease of maintenance for engineers and service technicians



## Features:

- » Generous weight capacity of 272kg to accommodate heavier users
- » Sleek and elegant yet robust design with secure lockable castors to maintain stability when in use
- » Electrically operated lifting and leg opening, reducing physical strain to caregivers
- » Excellent lift range from 50-170cm, providing effortless transfers, from floor to higher level surfaces
- » Large, braked rear castors for easy manoeuvring around furniture and over uneven flooring
- » Six-point cradle with spring clips offers a wider footprint, optimising user comfort and security during transfers
- » Ergonomic push handles offering a variety of comfortable hand positions
- » User-friendly handset
- » Battery can be charged on or off the hoist. Battery and external charger stand included as standard
- » LCD displays battery fuel indicator on patient handset and Lift Control Module
- » Emergency stop and manual lowering functions
- » Audible alerts to indicate low battery or weight overload
- » Maintenance monitoring: by months or lift cycles
- » Compatible with a complete range of Drive DeVilbiss slings



Six-point cradle



LCD information display



Electric leg opening



Effortless transfers from floor level



Large, braked rear castors



User-friendly handset



A complete range of Drive DeVilbiss slings are available to use with the Enhance 270 Hoist



Luxury Toilet Sling



Classic Universal Sling



In-Chair Hammock Comfort Sling

**Product Code:**

Enhance 270 Hoist	<b>PR51FLP600</b>
-------------------	-------------------

**Specifications:**

Safe Working Load	<b>272kg (600 lbs)</b>
Total Product Weight	<b>71kg (157 lbs)</b>
Castors	<b>7.6cm (Front) 12.5cm (Rear)</b>
Brakes	<b>Yes – Rear</b>
Base Length	<b>143.5cm</b>
External Leg Width (closed - open)	<b>77cm / 114cm</b>
Internal Leg Width (closed / open)	<b>62.5cm / 100cm</b>
Overall Height (min-max)	<b>68 - 194cm</b>
Lift Range	<b>50 - 170cm</b>
Turning Diameter	<b>142cm</b>
Height of Legs	<b>12cm</b>
Lifting/Lowering from control box	<b>Yes</b>

Overload Sensor	<b>Yes</b>
Motor Anti-Entrapment	<b>Yes</b>
Emergency Stop	<b>Yes</b>
Emergency Lowering/Lifting	<b>Yes - Manual/ Electric</b>
Battery	<b>Lead-Acid 2.9Ah</b>
Electrical Shock Protection	<b>TYPE BF - Class II</b>

For sling compatibility information please see our website [www.drivedevilbiss.co.uk](http://www.drivedevilbiss.co.uk)



info@drivedevilbiss.co.uk  
0845 0600 333

Leaflet Code: LL519

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.