

# LED Hydra

Unique, this new special "beauty care" lamp model allows the user to position the two lamps on either side of their hands, this avoids the production of shadows that occur when using only one light source.

Our Hydra lamp produces an scialitycal effect locally allowing minute work to be performed without shadows.

The lamp heads are mounted on two individually positionable flexible arms. A larger flexible arm connects the head unit to the base, providing support and manoeuvrability for ease of use.

Ideal in cosmetic surgery and beauty care for treating the face, scalp, for mesotherapy, hair implants, tattooing, etc.



HYDRA  
on/off Switch



● Medical device class 1

Light intensity **420 Lumens\* X 2**

Colour of temperature **4500°K\* (neutral white)**


Colour rendering index **CRI 95/100**

Electric consumption **6.5 Watts**

- LED life span: 50 000 h minimum
- Electrical protection: insulation class II
- Standards: IEC 60 601 -1/ -2-41
- Power cable length: 3 meters
- Suggested colors: white and black
- Head dimensions: 6 x 10 cm
- Flexible arm: 100 cm + 2 x 25 cm

• Fixings to choose from p23-26



Lux*	Beam 25°	Distance	Illuminated diameter
20000		30 cm	15 cm
7000		50 cm	25 cm
1700		100 cm	35 cm

VALUES FOR EACH HEAD

\* Values accurate within +/- 10 %

## CODES

On/off Switch	Contactless switch
LEDBI122	LEDBI122-ISC

