




# Veterinary scale "Easy" with fold function

8320

- ▶ Foldable veterinary scale for small animals
- ▶ Perfect for mobile use
- ▶ Low space requirement, easy to transport
- ▶ Hold and tare function
- ▶ Energy-saving auto-off function



## Features

---

- ▶ Patented folding mechanism with magnetic lock
- ▶ Low space requirement, easy to transport
- ▶ Easy to read LC display
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Switchable to kg/lb units
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 2 x 1.5 V AAA alkaline batteries (included in scope of delivery)

## Products

---

Order number Variation	Approvable	Division g	Maximum load kg
8320.01.001 Veterinary scale for small animals with fold function	No	10	15

---

## Accessories

---

Order number  
Variation

Description

5040.01.001



Carry bag black, with flexible shoulder strap and two carrying handles

---